

# Bedside cabinet with overbed table – bilateral functionality

Contemporary bedside cabinets that offer versatility and durability for today's patient environment.



B03



B03



B03



B04



B04



B04

## Features and benefits

- A varied selection of drawer and door configurations, suitable for all budgets
- Drawers that glide effortlessly into place
- Long-lasting, durable construction, using materials that are both easy to clean and maintain
- Aptly sized to suit the needs of today's patient environments



Optional: Plastic drawer insert



Optional: 2x single bottle holder for B03



Optional: 2x single bottle holder for B04



Optional: Towel holder (opposite to overbed table)



Optional: contrasting colours front (drawers and doors)



Optional: contrasting colours front and overbed table



Optional: contrasting colours front and overbed table



Optional: contrasting colours front (drawers and doors)

### Technical Specifications model SER - CH770A

Height incl. rails	870 mm
Width	605 mm
Depth	450 mm
Weight	~ 33 kg
<b>Overbed table</b>	
Height adjustment	610 – 950 mm
Tilting range	max. 60°
Surface area standard	547 x 349 mm
optional	617 x 349 mm
<b>Material</b>	
Body / Top plate	High Pressure Laminate (HPL), Melamine
Edges	Aluminium
Overbed table	High Pressure Laminate (HPL)

Castors	5 x 65 mm with 2 brakes
Bumper colour	grey
Handles/Colours	see handles and colours sheet
<b>Additional options</b>	
<ul style="list-style-type: none"> <li>• plastic drawer insert</li> <li>• towel holder</li> <li>• 2x single bottle holder</li> <li>• contrasting colours for                             <ul style="list-style-type: none"> <li>- body / top plate and /or overbed table</li> <li>- front (drawers and doors) and /or overbed table</li> <li>- overbed table</li> </ul> </li> </ul>	

Hill-Rom is a leading global medical technology company with more than 10,000 employees in over 100 countries. We partner with health care providers by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Clinical Workflow, Surgical Safety and Efficiency, and Respiratory Health.** Hill-Rom people, programs, and product brands work towards one mission: **Every day, around the world, we enhance outcomes for our patients and their caregivers.**

Not all products/options are available in all countries. For further information about this product or a service, please contact your local Hill-Rom representative or visit our webpage:

[www.hill-rom.com](http://www.hill-rom.com)

This document is destined solely for use by healthcare professionals. Devices shown in this brochure are intended for use with patients in departments of healthcare establishments. Hill-Rom recommends you carefully read the detailed instructions for safe and proper use included in the documents accompanying the devices. The personnel of the healthcare establishments are responsible for the proper use and maintenance of these devices.

Furniture manufacturer:  
Hill-Rom SAS, B.P. 14 - Z.I du Talhouët, 56330 Pluvigner - France (also European Representative)

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2016 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.

Doc. No: 5EN150320-03, 27 September 2016